BRASSELER

USA[®] DENTAL INSTRUMENTATION

High Speed Air Turbine Handpiece

NLXX

Rx Only

Please read this Operation Manual carefully before use and file for future reference.

1. User and Indications for Use

User: Qualified Professionals (Dentists)

Indications for Use: The NLXK Turbine is an air-powered dental handpiece intended for removing carious material, reducing hard tooth structure, cavity preparation, finishing tooth preparations and restorations and polishing teeth.

2. Precautions for handling and operation

- Please read these precautions carefully and use only as intended or instructed.
- Safety instructions are intended to avoid potential hazards that could result in personal injury or damage to the device. Safety instructions are classified as follows in accordance with the seriousness of the risk.

Class	Degree of Risk
 WARNING	Hazard that could result in serious injury or damage to the device if the safety instructions are not correctly followed.
 CAUTION	Hazard that could result in light or moderate injury or damage to the device if the safety instructions are not correctly followed.
NOTICE	General product specification information highlighted to avoid product malfunction and performance reduction.

WARNING

- Depressing the Push Button while the handpiece is in operation may cause overheating, serious technical damage and possible premature handpiece failure. During operation avoid contact with any oral tissue that may cause the Push Button to be depressed while the handpiece is in operation.

CAUTION

- Read this Operation Manual before use to fully understand the product functions and file for future reference.
- When operating the product always consider the safety of the patient.
- Users are responsible for the operational control, maintenance and continual inspection of this product.
- Do not attempt to disassemble the product nor tamper with the mechanism except as recommended by Brasseler USA in this Operation Manual.
- Do not allow any impact on to the product. Do not drop the product.
- Operators and all others in the area must wear eye protection and a mask when operating this handpiece.
- Should the product function abnormally, cease operation immediately and contact Brasseler USA.
- Do not use high acid water or sterilizing solutions to wipe, immerse or clean the product.

 **CAUTION**

- The following products are delivered in a non-sterile condition and must be autoclaved prior to use.
- Perform regular function and maintenance checks.
- If the product is not used for a long period check it is functioning correctly before using on a patient.
- To avoid clinical downtime it is recommended that a spare be kept on hand in case of a breakdown during treatment.
- U.S. Federal law restricts this device to sale by or on the order of a licensed physician.

3. Accessory List

No.	Part Name	Quantity	Remarks
1	Spray Nozzle	1	-
2	Cleaning Wire	1	Brush and wire set
3	PTL O-ring Set	1	Included in the package for NLXK-S and NLXK-T

4. Setting of Air & Water Supply Pressure

Measure the supply pressure at the handpiece/hose connection point and set the pressure to the value specified on the specification table (Fig. 1).

For Multi Gauge information refer to Option Parts List.


Fig. 1

WARNING

- Do not exceed the optimum pressure specified on the specification table.

CAUTION

- Do not use air contaminated by dust, moisture or oil.

5. Connection & Disconnection of the Handpiece

Refer to Operation Manuals of coupling and hose before connecting the handpiece.

NLXK-S and NLXK-T

5-1 Connection

- 1) Insert the handpiece into the Coupling (Fig. 2).
- 2) Make sure the handpiece is firmly connected to the coupling.

5-2 Disconnection

Pull back the Retention Lock Ring and remove the handpiece from the Coupling (Fig. 2).


Fig. 2

CAUTION

- Do not operate the Retention Lock Ring while under drive air pressure. The high pressure may cause sudden release of the handpiece from the coupling.

6. Insertion & Removal of the Bur

6-1 To Insert the Bur (Fig. 3)

- 1) Insert the bur into the chuck.
- 2) Depress the Push Button to open the chuck (①).
- 3) Fully insert the bur into the chuck until it stops (②) then release the Push Button.
- 4) Ensure that the bur is secure by gently pulling and pushing the bur WITHOUT depressing the Push Button.


Fig. 3

6-2 To Remove the Bur (Fig. 4)

- Depress the Push Button to open the chuck (①) and remove the bur (②).


Fig. 4

NOTICE

- Grip the handpiece while placing your thumb tip on the Push button which makes it easier to depress the button.

CAUTION

- Always insert the bur all the way into the chuck.
- Remove the bur only after the handpiece has completely stopped rotating.
- Always keep the bur shank clean.
- Entry of debris into the chuck, via the bur shank, could cause bur rotation slip and also prevent the bur from being securely located in the chuck.
- Do not exceed the bur speed recommended by the bur manufacturer.
- Do not exceed maximum bur length recommended by the handpiece manufacturer.
- Do not apply excess pressure to the bur as it may break or bend or become difficult to remove.
- DO NOT use burs with problems listed below as the bur may break, seize up or disengage from the chuck.
 - Bent, deformed, anisomorous (worn), rusted, broken, deficient bur.
 - Bur which is cracked on the edge or axis.
 - Non-ISO standard, or tampered bur.

7. Check before treatment

Check that the Head Cap is firmly tightened. Also check for handpiece vibration, noise and overheating outside the patient's oral cavity. If any abnormalities are found do not use the handpiece and contact Brasseler USA.

8. Maintenance

After each patient maintain the product as follows. Lack of maintenance could cause premature failure or overheating of the handpiece.

8-1 Cleaning handpieces with the Clean-Head System

After the treatment of each patient, clean the Clean Head.

- 1) Half fill a cup with clean water.
- 2) Remove dirt and debris from the Clean Head Holes with the cleaning wire and brush (supplied as an accessory for the handpiece) (Fig. 5).
- 3) Rotate the handpiece and immerse half of the handpiece head in the cup of water (Fig. 6).
- 4) Rotate then stop intermittently the handpiece 3 times for 2 to 3 seconds each time.
- 5) Wipe the handpiece dry.

* If the dirt could not be removed from the hole, clean it by brush.


Fig. 5


Fig. 6

8-2 Cleaning (Handpiece)

- 1) Remove debris from the product. Do NOT use a wire brush.
- 2) Wipe clean with alcohol-immersed cotton swab or cloth.

 This icon denotes that the product can be washed via Thermo Disinfector.

Refer to the Thermo-Disinfector manual.

CAUTION

- After washing with Thermo-Disinfector and prior to lubrication, dry the product until all internal moisture is thoroughly removed. Thermo-Disinfector moisture remaining inside the product could reduce the effect of lubrication and could cause corrosion inside of the product.
- To clean the product never use any solvent such as benzine or thinner.

8-3 Cleaning (Optic)

Wipe clean the Glass Rod tip with an alcohol-immersed cotton swab. Remove all debris and oil (Fig. 7).

CAUTION

- Do NOT use a sharp tool to clean the Glass Rod. It could damage the glass and reduce the light transmission.


Fig. 7

8-4 Cleaning (Spray Port)

When spray ports are clogged, or spray does not exit evenly from the ports, clean the ports (Fig. 8).

CAUTION

- Do not forcibly insert the wire into the port. Resultant port damage could cause the spray to be directed away from the bur, and the cooling efficiency reduces.
- Do not blow air into the Clean Head Hole.


Fig. 8

8-5 Lubrication

■ PANA SPRAY Plus or other approved spray lubricant

Apply PANA SPRAY Plus or other approved spray lubricant every time after each use and/or before autoclaving.

- 1) Remove the bur from the handpiece.
- 2) Insert the Spray Nozzle into the Spray Port nozzle on the can.
- 3) Insert the Spray Nozzle into the rear of the handpiece. Hold the handpiece and spray for approximately 2-3 seconds. Apply lubricant until it expels from the handpiece head for at least 2 seconds (Fig. 9).


Fig. 9

CAUTION

- When applying spray be sure to hold the handpiece firmly to prevent the handpiece from slipping out of the hand due to the spray pressure.
- Hold the spray can upright.

■ Chuck cleaning

Clean the Push Button chuck once a week.

- 1) Mount the Tip Nozzle into the spray can port.
- 2) Lubricate the chuck directly through the bur insertion hole (Fig. 10).
- 3) Lubricate handpiece by using PANA SPRAY Plus or other approved lubricant (Fig. 9) or Brasseler USA automatic handpiece cleaning and lubrication system.


Fig. 10

⚠ CAUTION

- If the chuck is not regularly cleaned the chuck grip may be weakened and the bur may be accidentally released while in operation.

■ Automatic handpiece cleaning and lubrication system

When using an automatic handpiece cleaning and lubrication system, refer to the system instructions.

8-6 Sterilization

Sterilize the product by autoclave sterilization. Remove the bur after each patient and sterilize as noted below.

- 1) Insert into an autoclave pouch. Seal the pouch.
- 2) Autoclavable under the conditions below.
Autoclave for more than 20 min. at 121°C, or 15 min. at 132°C, or 3 min. at 134°C.
- 3) The handpiece should remain in the autoclave pouch until required for use.

CAUTION

- Do not autoclave the product with other instruments even when it is in a pouch. This is to prevent possible discoloration and damage to the product from chemical residue on other instruments.
- Keep the product in suitable atmospheric pressure, temperature, humidity, ventilation, and sunlight. The air should be free from dust, salt and sulphur.
- Immediately after use, the product should be cleaned, lubricated and sterilized. If blood remains on the external or internal surfaces it can become clotted and cause rust.
- Do not heat or cool the product too quickly. Rapid change in temperature could cause damage to the product.
- If the sterilizer chamber temperature may exceed 135°C during the dry cycle then delete the dry cycle.

 **CAUTION**

- Autoclave sterilization is recommended for the product. The validity of other sterilization methods is not confirmed.
- Do not touch the product immediately after autoclaving as it will be very hot and must remain in a sterile condition.

NOTICE

- Brasseler USA recommends Class B sterilizers as stated in EN13060.

9. Replacing the O-rings (NLXK-S and NLXK-T)

Replace the O-rings if water is present in the exhaust air line. This is an indication of possible water leakage within the coupling. ALWAYS change the complete set of O-rings.

- 1) Loosen and remove the taper ring at the rear of the handpiece (Fig. 11).
- 2) Gently remove each O-ring by hand (Fig. 12).
- 3) Insert the complete set of new O-rings in the correct grooves.
- 4) Firmly tighten the taper ring.

* Refer to Spare Parts List to identify the correct parts.

CAUTION

- Do not force the new replacement O-ring with excessive pressure.
- When inserting new O-rings, make sure they are inserted in the correct grooves.
- Make certain that the taper ring is firmly tightened. If the taper ring is loose water and air leakage could occur.


Fig. 11


Fig. 12

10. Replacing the Non-Retraction Valve

A water Non-Retraction Valve is integrated in the Coupling Joint, which shuts off the water retraction directly at the handpiece head to prevent fluids infiltrating the water line. If water is beginning to leak from the handpiece, replace the non-retraction valve.

■ Brasseler USA Coupling (NLXK-S and NLXK-T)

- 1) Remove the Coupling Joint from the hose.
 - 2) Remove the back-end Gasket.
 - 3) Pull and remove the water tube, and remove the Non-Retraction Valve.
 - 4) Insert the new Non-Retraction Valve securely and remount the back-end Gasket.
- * Refer to Spare Parts List to identify the correct parts.


Fig. 13

11. Periodical Maintenance Checks

Perform periodical maintenance checks every three months, referring to the check sheet below. If any abnormalities are found, contact Brasseler USA.

Points to check	Details
Head cap is loose	Check that the Head Cap is firmly tightened.
Rotation	Rotate the handpiece and check for abnormalities such as abnormal rotation, vibration, noise, and overheating.
Coolant Water	Operate the handpiece and check that the coolant water is flowing through all spray ports.

12. Symbol

This product is Autoclavable up to Max.135°C. (ISO 21531)

This product can be washed via Thermo Disinfector. (ISO 21531)

Manufacturer. (ISO 15223-1)

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed physician.

(FDA 21 CFR1 Part 801.109 (b)(1))

Catalogue Number. (ISO 15223-1)

QTY Quantity.

13. Warranty

Brasseler USA products are warranted against manufacturing errors and defects in materials. Brasseler USA reserves the right to analyze and determine the cause of any problem. Warranty is voided should the product be not used correctly or for the intended purpose or has been tampered with by unqualified personnel or has had non Brasseler USA parts installed. Replacement parts are available for seven years beyond discontinuation of the model.

14. Option Parts List

Model	Order Code
MG-4H Multi Gauge	Z109400

15. Spare Parts List

Model	Order Code	Compatible Product	Remarks
PTL O-ring Set	Y900580	NLXK-S and NLXK-T	-
Non-Retraction Valve	P401054	Brasseler USA coupling	-
Cleaning Wire	Z144095	-	-

* For other parts contact Brasseler USA.

16. Disposing product

In order to avoid the health risks of operators handling the disposal of medical equipment, as well as the risks of environmental contamination caused thereof, a surgeon or a dentist is required to confirm the equipment is sterile. Ask specialist firms who are licensed to dispose of specially controlled industrial wastes, to dispose the product for you.

17. Specifications

Model	NLXK-S	NLXK-KS
Rotation Speed	360,000 - 440,000min ⁻¹ (at the recommended drive air pressure)	
Bur Type	ISO 1797-1 Ø1.59 - 1.60mm Short Shank Bur/Standard Bur	
Chucking Length	10mm	
Max. Bur Length	21mm	
Max. Working Part Diameter	Ø2mm	
Recommended Drive Air Pressure	0.28MPa (2.8kgf/cm ²)	
Operating Drive Air Pressure Range	0.25 - 0.30MPa (2.5 - 3.0kgf/cm ²)	
Max. Air Consumption	<66NL/min (0.3MPa)	
Water Consumption	≥45mL/min (0.1MPa)	
Chip Air Consumption	≥1.5NL/min (0.15MPa)	
Water Pressure	0.08 - 0.20MPa (0.8 - 2.0kgf/cm ²)	
Chip Air Pressure	0.10 - 0.20MPa (1.0 - 2.0kgf/cm ²)	
Optic	Glass Rod	
Use Environment	Temperature: 10 - 40°C (No Condensation), Humidity: 30 - 75%	
Transportation and Store Environment	Temperature: -10 - 50°C, Humidity: 10 - 85%, Atmospheric Pressure: 500 - 1,060hPa	

Model	NLXK-T	NLXK-KT
Rotation Speed	320,000 - 400,000min ⁻¹ (at the recommended drive air pressure)	
Bur Type	ISO 1797-1 Ø1.59 - 1.60mm Standard Bur	
Chucking Length	12mm	
Max. Bur Length	21mm	
Max. Working Part Diameter	Ø2mm	
Recommended Drive Air Pressure	0.28MPa (2.8kgf/cm ²)	
Operating Drive Air Pressure Range	0.25 - 0.30MPa (2.5 - 3.0kgf/cm ²)	
Max. Air Consumption	<66NL/min (0.3MPa)	
Water Consumption	≥45mL/min (0.1MPa)	
Chip Air Consumption	≥1.5NL/min (0.15MPa)	
Water Pressure	0.08 - 0.20MPa (0.8 - 2.0kgf/cm ²)	
Chip Air Pressure	0.10 - 0.20MPa (1.0 - 2.0kgf/cm ²)	
Optic	Glass Rod	
Use Environment	Temperature: 10 - 40°C (No Condensation), Humidity: 30 - 75%	
Transportation and Store Environment	Temperature: -10 - 50°C, Humidity: 10 - 85%, Atmospheric Pressure: 500 - 1,060hPa	

* Available models differ according to countries/regions.

NLXK-S / NLXK-T

For Brasseler USA Phatelus/Mach Coupling

NLXK-KS / NLXK-KT

For KaVo® MULTiflex®, MULTiflex® LUX Coupling

KaVo® and MULTiflex® are registered trademarks of Kaltenbach & Voigt GmbH & Co., Germany.

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KaVo® et MULTiflex® sont des marques déposées de la société Kaltenbach & Voigt GmbH & Co., Allemagne.
Brasseler U.S.A. Dental, LLC
One Brasseler Blvd Savannah, GA 31419 U.S.A.