**BRASSELER**  
USA®

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**FORZA+ X45**

## 1 User and Indications for Use

User: Licensed dentists / oral surgeons

Indications for Use: FORZA+ X45 is powered by either an air-motor or electronic micromotor for use in general dentistry. The device is intended for cutting and grinding teeth, cavity preparations, tooth and crown preparations, finishing and trimming teeth and filling materials, and removal of crowns and filling materials.

## 2 Precautions for handling and operation

- Please read these precautions carefully and use only as intended or instructed.
- Safety instructions are intended to avoid potential hazards that could result in personal injury or damage to the device. Safety instructions are classified as follows in accordance with the seriousness of the risk.

Class	Degree of Risk
 <b>WARNING</b>	Hazard that could result in serious injury or damage to the device if the safety instructions are not correctly followed.

Class	Degree of Risk
<b>⚠ CAUTION</b>	Hazard that could result in light or moderate injury or damage to the device if the safety instructions are not correctly followed.
<b>NOTICE</b>	General product specification information highlighted to avoid product malfunction and performance reduction.

## **⚠ WARNING**

- Immediately after a treatment, the handpiece should be cleaned, lubricated and sterilized. Failure to properly maintain the handpiece may lead to overheating, causing burn injuries or product failure. (Refer to “8. Post-use Maintenance”)
- If blood infiltrates inside a handpiece, an automatic handpiece cleaning and lubrication system may not totally clean the internal handpiece components, which may lead to internal blood coagulation. Coagulated blood inside the handpiece may cause handpiece failure and overheating, causing burn injuries. To avoid the risk, Brasseler USA recommends to immediately lubricate those handpieces with PANA SPRAY Plus.

**⚠ WARNING**

- Depressing the Push Button while the handpiece is in rotation may lead to overheating, causing burn injuries or product failure. Avoid the push button to contact with any oral tissue.
- Supply coolant water and coolant air while using the handpiece. No supplying the coolant water and coolant air may lead to overheating, causing burn injuries or product failure.
- Do not use the handpiece for polishing purposes. Polishing paste could enter the handpiece, causing the failure of the Push Button or burn injuries due to the overheating of the handpiece head.

**⚠ WARNING**

- Keep any debris or other foreign materials away from inside the gear or the handpiece. Foreign materials remaining inside may lead to overheating, causing burn injuries or other accidents.

## **WARNING**

- Ball bearings are wear items and may require replacement. Be sure to inspect the handpiece before use (Refer to “7. Checking the Handpiece Before Each Use”). If abnormal vibration or noise are found, the bearings or other internal parts may be worn out or damaged. Damaged parts may lead to handpiece overheating causing burn injuries. If overheating is suspected, stop using the handpiece immediately and contact Brasseler USA.

## **CAUTION**

- Read this Operation Manual before use to fully understand the product functions and file for future reference.
- When operating the product always consider the safety of the patient.
- Users are responsible for the operational control, maintenance and continual inspection of this product.
- Prior to clinical use, inspect the handpiece. Check for vibration, noise and overheating. If any abnormalities are found, stop using the handpiece immediately and contact Brasseler USA. (Refer to “7. Checking the Handpiece Before Each Use”)

## ⚠ CAUTION

- Do not disassemble or alter the handpiece except as recommended by Brasseler USA in this Operation Manual.
- Do not allow any impact on to the product. Do not drop the product.
- Operators and all others in the area must wear eye protection and a mask when operating this handpiece.
- Do not use burs with problems listed below as the bur may break, seize up or disengage from the chuck.
  - Bent, deformed, worn, rusted, broken, deficient bur.
  - Bur which is cracked on the edge or axis.
  - Non-ISO standard, or tampered bur.

## ⚠ CAUTION

- Do not use burs other than specified in “10. Specifications” as such burs may be accidentally released or break during rotation.
- Always keep the bur shank clean. Dirt or debris in the chuck could cause poor bur concentricity or low chuck retention force.
- Always insert the bur all the way into the chuck. If insertion is insufficient, premature failure of the bearings or accidental release of the bur may occur.
- Do not exceed Max. Bur Length shown in “10. Specifications”.

## **⚠ CAUTION**

- Do not use a surgical bur with its maximum length of 25mm or longer. Using such bur with this handpiece may lead to premature abrasion of the bearings. In addition, sudden release, bend or breakage of the bur may occur.
- Always follow the instructions provided by the bur manufacturer.
- Do not exceed the bur speed recommended by the bur manufacturer.
- Should the handpiece function abnormally during use, stop using the handpiece immediately and contact Brasseler USA.

## **⚠ CAUTION**

- Do not use the following fluids to wipe, immerse or clean the product; strong/super acid water, strong acid/alkaline chemicals, chlorine-containing solutions, solvents such as benzene or thinner.
- The handpiece is delivered in a non-sterile condition and must be autoclaved prior to use.
- Perform regular function and maintenance checks. (Refer to “9. Periodical Maintenance Checks”)
- If the handpiece has not been used for a long period, rotate the handpiece and check for noise, vibration and overheating before use.
- To avoid clinical downtime it is recommended that a spare be kept on hand in case of a breakdown during treatment.

**⚠ CAUTION**

- The operation of the handpiece is permitted only on dental units which correspond to the standards IEC 60601-1 (EN 60601-1) and IEC 60601-1-2 (EN 60601-1-2).
- U.S. Federal law restricts this device to sale by or on the order of a licensed physician.

**3 Accessory List**

No.	Part Name	Quantity
1	Wrench	1
2	E-Type Spray Nozzle	1

## 4 Connecting and Disconnecting the Handpiece and the Motor

### 4-1 Connecting

- 1 Insert the handpiece direct to the motor (Fig. 1). Screw the handpiece until it clicks and locks into position.
- 2 Confirm that the handpiece is firmly connected to the motor.


Fig. 1

### 4-2 Disconnecting

Hold the motor and the handpiece, then pull apart.

#### **⚠ CAUTION**

- Do not connect or disconnect the handpiece until the motor has completely stopped.
- Connect ONLY to E type motors (ISO 3964).
- Do not exceed the Max. Rotation Speed (Motor) shown in "10. Specifications".

## 5 Mounting and Removing the Bur

### 5-1 Mounting the Bur

- 1 Insert the bur into the chuck.
- 2 Depress the Push Button to open the chuck (1).
- 3 Insert the bur fully into the chuck until it stops (2) then release the Push Button.
- 4 Confirm that the bur is not removed by gently pulling and pushing the bur without depressing the Push Button.


Fig. 2

### NOTICE

- Grip the handpiece while placing your thumb tip on the Push button which makes it easier to depress the button.

## 5-2 Removing the Bur

Depress the Push Button to open the chuck (①) and remove the bur (②).


Fig. 3

### ⚠ CAUTION

- Do not use burs with problems listed below as the bur may break, seize up or disengage from the chuck.
  - Bent, deformed, worn, rusted, broken, deficient bur.
  - Bur which is cracked on the edge or axis.
  - Non-ISO standard, or tampered bur.
- Do not use burs other than specified in "10. Specifications" as such burs may be accidentally released or break during rotation.
- Always keep the bur shank clean. Dirt or debris in the chuck could cause poor bur concentricity or low chuck retention force.

**⚠ CAUTION**


- Always insert the bur all the way into the chuck. If insertion is insufficient, premature failure of the bearings or accidental release of the bur may occur.
- Do not exceed Max. Bur Length shown in “10. Specifications”.
- Do not use a surgical bur with its maximum length of 25mm or longer. Using such bur with this handpiece may lead to premature abrasion of the bearings. In addition, sudden release, bend or breakage of the bur may occur.
- Always follow the instructions provided by the bur manufacturer.

**⚠ CAUTION**

- Do not exceed the bur speed recommended by the bur manufacturer.
- Do not mount or remove the bur until the motor has completely stopped.
- Do not apply excess pressure to the bur as it may break or bend or become difficult to remove.
- When using a bur with a length of 20-25mm, the working diameter must be within the value specified in the table below. The use of a bur with a diameter that exceeds the specified values may result in a bur run out and cause the bur to bend or break.

## ⚠ CAUTION

øa	Diamond Bur: ø2mm or less
	Carbide Bur : ø1mm or less

Be careful of overheating when using some types of carbide bur, as the coolant water (in "water jet" mode) cannot reach the bur tip because of the shape of the tip.

## ⚠ CAUTION

Brasseler USA recommends carbide burs shown in the table below.

Manufacturer	Shape	No.	Max. Rotation Speed (min <sup>-1</sup> )
Midwest	ROUND	#1	300,000
	STRAIGHT FLAT ENDCROSS CUT	#556	300,000
Maillefer	CYLINDER SQUARE	#557L	300,000
MANI	ROUND	#2 (L=25)	120,000
	PEAR	#330 (L=25)	120,000
	CYLINDER ROUND	#1557 (L=25)	120,000
BRASSELER	FG SURG FISSURE	#557 (L=25)	300,000

\*Adjust the handpiece's speed so that the speed of the bur does not exceed the specified maximum value.

## 6 Switching coolant water type

It is possible to switch the type of coolant water flow which cools the drilling point.

Turn the switching valve with an attached wrench to switch the coolant water type.

- Align the mark  and  to set the water spray shape.
- Align the mark  and  to set the water jet shape.


Fig. 4

## CAUTION

- Make sure to turn the switching valve until it stops.
- Before initiating treatment, assure the water is spraying correctly.
- If water is not being supplied correctly, recheck the switching valve.

## 7 Checking the Handpiece Before Each Use

---

Follow the check below before use. If any abnormalities are found, stop using the handpiece immediately and contact Brasseler USA.

- ➊ Check the Head Cap is firmly tightened.
- ➋ Check the coolant water is flowing properly.
- ➌ Mount the bur. (Refer to “5. Mounting and Removing the Bur”)
- ➍ Rotate the handpiece for about one minute with coolant water at the Max. rotation speed of the attached motor. During rotation, check for abnormalities such as abnormal rotation, vibration, noise.

- 5 After the handpiece rotation has completely stopped, touch the handpiece head to confirm the head is NOT heating abnormally (Fig. 5).


Fig. 5

**⚠ CAUTION**

- To avoid injury, keep your hands away from the bur during rotation.

## 8 Post-use Maintenance

After each patient maintain the product as follows.

**⚠ WARNING**

- Immediately after a treatment, the handpiece should be cleaned, lubricated and sterilized. Failure to properly maintain the handpiece may lead to overheating, causing burn injuries or product failure.

**⚠ CAUTION**

- Do not use the following fluids to wipe, immerse or clean the product; strong/super acid water, strong acid/alkaline chemicals, chlorine-containing solutions, solvents such as benzine or thinner.

## 8-1 Preparation

- 1 Wear a eye protection, a mask and gloves to avoid infection.
- 2 Remove the bur.
- 3 Disconnect the handpiece from the motor.

## 8-2 Cleaning

### ■ Manual Cleaning (External Surface)

- 1 Remove dirt and debris on the handpiece surface using a brush under running water (at 38°C (100°F) or less). Do not use a wire brush.
- 2 Wipe off the moisture on the surface using absorbent cloth.

### CAUTION

- If water comes inside the handpiece, blow air into the rear of the handpiece using an air syringe and remove the remaining water.

■ Automatic Cleaning (External and Internal Surfaces)  
Use a thermo-disinfector in accordance with ISO15883-1.  
Refer to the thermo-disinfector manual.

**⚠ CAUTION**

- After washing with thermo-disinfector and prior to lubrication, dry the product until all internal moisture is thoroughly removed. Thermo-disinfector moisture remaining inside the product could reduce the effect of lubrication and could cause corrosion inside of the product.

**8-3 Cleaning (Optic)**

Wipe clean the Glass Rod tip with an alcohol-immersed cotton swab. Remove all debris and oil.


Fig. 6

**⚠ CAUTION**

- Do not use a sharp tool to clean the Glass Rod. It could damage the glass and reduce the light transmission.

## 8-4 Disinfection

### ■ Manual Disinfection (External Surface)

Wipe the surface with cloth soaked in disinfectant solution such as alcohol.

### ■ Automatic Disinfection (External and Internal Surfaces)

Use a thermo-disinfector in accordance with ISO15883-1. Refer to the thermo-disinfector manual.

## 8-5 Lubrication

### ■ PANA SPRAY Plus

Apply PANA SPRAY Plus every time after each use and/or before autoclaving.

- 1 Attach the E-Type Spray Nozzle into the Spray Port nozzle on the can.

- 2 Insert the E-Type Spray Nozzle into the rear of the handpiece. Hold the handpiece and spray for approximately 2-3 seconds. Apply lubricant until it expels from the handpiece head for at least 2 seconds (Fig. 7).


Fig. 7

## ⚠ CAUTION

- Hold the spray can upright.
- Be sure to hold the handpiece firmly to prevent slipping out of hand by the spray pressure when lubricating.
- Supply lubricant until it comes out of the handpiece head.
- If you wish to expel excessive oil from inside of the handpiece, rotate the handpiece for approx.15 seconds without a bur. During the rotation, do not depress the push button inadvertently, especially when using absorbent cloth to prevent oil from scattering. Depressing the push button during rotation may cause a failure to attach a bur.

## ⚠ CAUTION


Fig. 8

Fig. 9

## NOTICE

- Brasseler USA recommends the use of “Spray Mist Absorber” (REF Y900084) to prevent oil mist expelling out of the handpiece head.

### ■ Chuck cleaning

Clean the Push Button chuck once a week.

- 1 Mount the Tip Nozzle into the spray can port.
- 2 Lubricate the chuck directly through the bur insertion hole (Fig. 10).
- 3 Lubricate handpiece by using PANA SPRAY Plus (Fig. 7) or automatic handpiece cleaning and lubrication system.


Fig. 10

### ⚠ CAUTION

- If the chuck is not regularly cleaned the chuck grip may be weakened and the bur may be accidentally released during use.

■ Automatic handpiece cleaning and lubrication system  
When using automatic handpiece cleaning and lubrication system, refer to the system instructions.

## **WARNING**

- If blood infiltrates inside a handpiece, an automatic handpiece cleaning and lubrication system may not totally clean the internal handpiece components, which may lead to internal blood coagulation. Coagulated blood inside the handpiece may cause handpiece failure and overheating, causing burn injuries. To avoid the risk, Brassler USA recommends to immediately lubricate those handpieces with PANA SPRAY Plus.

## **8-6 Sterilization**

Sterilize the product by autoclave sterilization. After each patient sterilize as noted below.

- ➊ Insert into an autoclave pouch. Seal the pouch.
- ➋ Autoclavable under the conditions below.
  - Gravity displacement steam sterilizer; Exposure at 132 °C (270 °F) for 15 minutes, Drying for 30 minutes/  
Exposure at 135 °C (275 °F) for 10 minutes, Drying for 30 minutes
  - Dynamic air removal steam sterilizer: Exposure at 132 °C (270 °F) for 4 minutes, Drying for 30 minutes/  
Exposure at 135 °C (275 °F) for 3 minutes, Drying for 16minutes
- ➌ The handpiece should remain in the autoclave pouch until required for use.

### **⚠ CAUTION**

- Be sure to clean and lubricate the handpiece prior to autoclaving. If blood remains on the internal surface it can become clotted and cause product failure.
- Do not use sterilization pouches containing water-soluble adhesive components such as PVA (polyvinyl alcohol). Adhesive components washed out during sterilization may enter the product, causing problems such as poor rotation, sticking and failure to operate correctly.
- Do not autoclave the product with other instruments even when it is in a pouch. This is to prevent possible discoloration and damage to the product from chemical residue on other instruments.
- Do not heat or cool the product too quickly. Rapid change in temperature could cause damage to the product.

### **⚠ CAUTION**

- To avoid product failure, do not use a sterilizer whose temperature exceeds 138°C during use including the dry cycle. The chamber temperature of some type of sterilizers may exceed 138°C. Contact the sterilizer manufacturer for detailed information.
- Keep the product in suitable atmospheric pressure, temperature, humidity, ventilation, and sunlight. The air should be free from dust, salt and sulphur.
- Do not touch the product immediately after autoclaving as it will be very hot and must remain in a sterile condition.
- Autoclave sterilization is recommended for the product. The validity of other sterilization methods is not confirmed.

## NOTICE

- Brasseler USA recommends Class B sterilizers as stated in EN 13060.

## 9 Periodical Maintenance Checks

Perform periodical maintenance checks every three months, referring to the check sheet below. If any abnormalities are found, contact Brasseler USA.

Points to check	Details
Head cap is loose	Check that the Head Cap is firmly tightened.
Rotation	Rotate the handpiece and check for abnormalities such as abnormal rotation, vibration, noise, and overheating.
Coolant Water	Operate the handpiece and check that the coolant water is flowing through all spray ports.

## 10 Specifications

Model	FORZA+ X45
Max. Rotation Speed (Motor)	40,000min <sup>-1</sup>
Max. Rotation Speed (Handpiece)	168,000min <sup>-1</sup>
Gear Ratio	1:4.2 Increasing
Bur Type	ISO 1797 Type3 1.59-1.60mm CA Bur
Chuck Type	Push Button Type
Chuck Length	11.6mm
Max. Bur Length	25mm
Max. Working Part Diameter	Ø2mm
Optic	Glass Rod

Coolant Water Type	Water Spray or Water Jet (switch-selectable)
Water Flow Rate	Min. 50ml/min (when the water pressure is 0.2MPa )
Chip Air Consumption	Min. 1.5NI/min (when the water pressure is 0.2MPa )
Use Environment	Temperature: 10 - 35°C (50- 95°F), Humidity: 30 - 75% (No Condensation)
Transportation and Store Environment	Temperature -10 - 50°C (14- 122°F), Humidity: 10 - 85%, Atmospheric Pressure: 500 - 1,060hPa

## 11 Symbol

---

The instrument can be used with a washer-disinfecter for thermal disinfection.

This product can be sterilized in a steam sterilizer at 135°C (275°F).

Consult operation instructions

Manufacturer

Medical Device

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed physician.

Device Name for Health Canada

Catalog number (Order code)

Serial No.

Date of manufacture

Unique Device Identifier

GS1 DataMatrix for Unique Device Identifier.

Temperature limitation

Humidity limitation

Atmospheric pressure limitation

With the function of the switching coolant water type.

## 12 Warranty

Brasseler USA products are warranted against manufacturing errors and defects in materials. Brasseler USA reserves the right to analyze and determine the cause of any problem. Warranty is voided should the product be not used correctly or for the intended purpose or has been tampered with by unqualified personnel or has had non Brasseler USA parts installed. Replacement parts are available for seven years beyond discontinuation of the model.

## 13 Spare Parts List

Model	Order Code	Remarks
Wrench	5025747U0	For switching coolant water type
E-Type Spray Nozzle	Z019090	-

## **14** Disposing product

---

In order to avoid the health risks of operators handling the disposal of medical equipment, as well as the risks of environmental contamination caused thereof, a surgeon or a dentist is required to confirm the equipment is sterile. Ask specialist firms who are licensed to dispose of specially controlled industrial wastes, to dispose the product for you.

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**Brasseler U.S.A. Dental, LLC**
One Brasseler Blvd, Savannah, GA 31419 U.S.A.