



BRASSELER USA®
DENTAL INSTRUMENTATION

One Brasseler Boulevard
Savannah, GA 31419
t: 800.841.4522
f: 888.610.1937

MEDIA CONTACT:

Erin Foster/Lanmark360
efoster@lanmark360.com
Tel: 732-389-4500 x128

FOR IMMEDIATE RELEASE:

Brasseler USA® Introduces HygienePro™ Ultrasonic Scaler Inserts

***24 Ergonomic Scaler Inserts Available for Greater Patient Comfort
and Practice Efficiency***

SAVANNAH, GA (September 10, 2014) – Brasseler USA®, a leading manufacturer of quality instrumentation, is pleased to introduce a wide variety of new ultrasonic scaler inserts for its popular HygienePro™ family of products for the dental hygiene professional.

HygienePro™ Ultrasonic Scaler Inserts are ergonomically designed to substantially reduce user fatigue while increasing patient comfort. Available in both hard and soft grips, the new inserts are extremely versatile and easy to implement. HygienePro Ultrasonic Scaler Inserts are compatible with all magnetostrictive devices.

Ergonomic grip designs feature a textured, non-slip surface that significantly reduces hand and wrist fatigue for enhanced user comfort and control. Moreover, HygienePro Ultrasonic Scaler Inserts are equipped with a patented laminate solder joint ensuring extraordinary durability, consistent spray, and an even vibration.

A wide variety of tip designs are available in both 25K and 30K frequency options to suit each user's unique individual requirements:

- **Universal** – the universal tip design is available in slotted and direct focus water delivery options
- **Slim** – the slim tip design is also available in slotted and direct focus water delivery options
- **Triple Bend** – the triple bend tip design is available exclusively in slotted water delivery option
- **Beaver Tail** – the beaver tail tip design is available exclusively in slotted water delivery option



BRASSELER USA®
DENTAL INSTRUMENTATION

One Brasseler Boulevard
Savannah, GA 31419
t: 800.841.4522
f: 888.610.1937

For more information about the HygienePro Ultrasonic Scaler Inserts, please visit www.BrasselerUSA.com or call 800-841-4522.

About Brasseler USA®

Brasseler USA® is a leading ISO Certified healthcare company, providing quality instrumentation to healthcare professionals for use in restorative dentistry, endodontics, oral surgery and oral hygiene. Over the past 35 years, Brasseler USA has developed a reputation as an innovative market leader in diamonds, carbides, polishers, endodontics, hand instruments and handpieces. Today, Brasseler USA offers the most comprehensive assortment of instruments and power systems under one brand in the world. For more information, please visit www.BrasselerUSA.com.

###