






















**BRASSELER**  
USA®

**NEW**  
**K0291**

Restorative Design Prep Kit  
DEVELOPED BY  
**DR. ROBERT WINTER**



										
<b>Winter No.</b>	<b>WS8M.6</b>	<b>WS8F.6</b>	<b>WS8M1.0</b>	<b>WS8F1.0</b>	<b>WS8M1.2</b>	<b>WS8F1.2</b>	<b>WS6M1.2</b>	<b>WS6F1.2</b>	–	<b>WSMR1.0</b>
<b>Figure No.</b>	<b>856</b>	<b>8856</b>	–	–	–	–	–	–	<b>8856</b>	–
Head Size (1/10 mm)	012	012	018	018	020	020	018	018	016	012
Cutting Length (mm)	8.0	8.0	8.0	8.0	8.0	8.0	6.0	6.0	8.0	0.27
Grit	Medium	Fine	Medium	Fine	Medium	Fine	Medium	Fine	Fine	Fine
Shank Type	FG	FG	FG	FG	FG	FG	FG	FG	FG	FG
Quantity	1	1	1	1	1	1	1	1	1	1

									
<b>Winter No.</b>	<b>RWMIN</b>	<b>RWCONV</b>	<b>RWEXT</b>	<b>WDL1.0</b>	<b>WDL1.5</b>	<b>WDL2.0</b>	<b>WDL2.5</b>	–	–
<b>Figure No.</b>	<b>868C</b>	<b>868D</b>	<b>868E</b>	–	–	–	–	<b>379</b>	<b>379F</b>
Head Size (1/10 mm)	026	030	034	016	017	017	018	023	023
Cutting Length (mm)	–	–	–	–	–	–	–	4.2	4.2
Depth Cut (mm)	.3 / .5 / .7	.5 / .7 / .9	.8 / 1 / 1.2	1.0	1.5	2.0	2.5	–	–
Grit	Medium	Medium	Medium	Medium	Medium	Medium	Medium	Medium	Fine
Shank Type	FG	FG	FG	FG	FG	FG	FG	FG	FG
Quantity	1	1	1	1	1	1	1	1	1

**BY YOUR SIDE**