



## SPLINT SYSTEM | Developed by Brasseler USA

Designed to adjust and polish splint appliances, this system consists of E-CUTTER® instruments for trimming, diamond discs to finish and contour, a Scotch-Brite™ buff to remove tags left behind by trimming and various acrylic polishers, including ET® ProviPro™ cotton buffs to restore surface smoothness. This complete system is the ideal solution when fitting your patients with splints.

### K0221

Figure No.	0674	0654	H251EF	H351EF	H79EF	365
Head Size (1/10 mm)	100	100	060	060	040	220
Cutting Length (mm)	24.0	24.0	14.0	11.0	13.0	-
Grit	Coarse	Fine	Fine	Fine	Fine	Medium
Shank Type	HP	HP	HP	HP	HP	HP
Quantity	1	1	1	1	1	1
Description	AcryliPro™ Polishing Bullet: Extra-oral polishing of acrylics		EF Cutter: Fine E-CUTTER® for trimming and finishing acrylic based appliances			Diamond Disc: Serrated for acrylics

Figure No.	SWR22M	SWR22F	RWGPP	25125HP	305
Head Size (1/10 mm)	220	220	220	250	050
Cutting Length (mm)	3.0	3.0	3.0	-	-
Grit	Medium	Fine	Medium	Fine	-
Shank Type	-	-	-	HP	HP
Quantity	1	1	1	1	3
Description	ET® ProviPro™ Polishing Buff: Polishing of acrylic and composite based appliances		Silicone Knife-Edge Polisher: Shaping and smoothing prior to final polish	Scotch-Brite™ Buff: Removal of tags created from cutters	Mounting for ET® ProviPro™ and Silicone Knife-Edge Polisher

Scotch-Brite™ is a trademark of 3M.

By Your Side  
in Dentistry



**BRASSLER USA**  
DENTAL INSTRUMENTATION



To order call 800.841.4522 or fax 888.610.1937.  
In Canada call 800.363.3838 or fax 866.330.4454.  
Visit our website at BrasselerUSA.com

