


solutions™

Dr. Ronald E. Goldstein
ET® Carbide Combo Kit



Order Number:
4200



US Number	ET3	ET3F	ET3UF	ET4	ET4F	ET4UF	ET6	ET6F	ET6UF	ET9
Figure Number	H132	H132F	H132UF	H133	H133F	H133UF	H134	H134F	H134UF	H135
Head Size ($\frac{1}{10}$ mm)	008	008	008	010	010	010	014	014	014	014
Head Length (mm)	3.0	3.0	3.0	4.0	4.0	4.0	6.0	6.0	6.0	9.0



US Number	ET9F	ET9UF	OS1	OS1F	OS1UF	OS2	OS2F	OS2UF	OS3	OS3F	OS3UF
Figure Number	H135F	H135UF	H379	H379F	H379UF	H247	H247F	H247UF	H247	H247F	H247UF
Head Size ($\frac{1}{10}$ mm)	014	014	023	023	023	009	009	009	007	007	007
Head Length (mm)	9.0	9.0	4.2	4.2	4.2	3.2	3.2	3.2	3.2	3.2	3.2

Only Available Direct from your Brasseler Representative

To order call 800.841.4522 or fax 888.610.1937

In Canada call 800.363.3838 or fax 866.330.4454













Visit our website at BrasselerUSA.com










By Your Side In Dentistry



B-3732B-REVI-5022192L0-01.13-RRD

solutions™ 4200 ET® Dr. Ronald E. Goldstein Carbide Combo Kit

	Order #	Head Size	Grit	Usage
	ET3	008	F	Large amount gingival/interproximal finishing of direct bonded restorations
	ET3F	008	EF	Minimal amount of gingival or proximal finishing for direct resin restorations
	ET3UF	008	UF	Final subgingival/interproximal finishing of bonded/provisional/ceramic restorations
	ET4	010	F	Gross gingival finishing of direct bonded restorations
	ET4F	010	EF	Minimal gingival finishing of direct bonded or provisional restorations is required
	ET4UF	010	UF	Final finishing of direct bonded, provisional, or ceramic restorations
	ET6	014	F	Gross labial/lingual finishing of medium sized teeth
	ET6F	014	EF	Labial/lingual minimal finishing anterior/posterior teeth
	ET6UF	014	UF	Final finishing or placing texture in labial surfaces of bonded or provisional restorations
	ET9	014	F	Gross contouring of labial surfaces of direct resin or provisional restorations
	ET9F	014	EF	Minimal contouring or finishing of labial or lingual surfaces on longer teeth
	ET9UF	014	UF	Texturing/final finish of labial surfaces on longer bonded/provisional restorations

	Order #	Head Size	Grit	Usage
	OS1	023	F	Gross lingual/occlusal finishing and occlusal equilibration procedures
	OS1F	023	EF	Minimal finishing lingual/occlusal surfaces of bonded/provisional restorations or slight adjustment of natural teeth
	OS1UF	023	UF	Final occlusal/lingual finishing bonded restorations/occlusal adjustments temporary provisional materials
	OS2	009	F	Occlusal anatomy in posterior bonded or provisional restorations
	OS2F	009	EF	Initial placement of occlusal anatomy when lesser amounts are required
	OS2UF	009	UF	Final polishing/finishing in posterior bonded/temporary restorations
	OS3	007	F	supplemental occlusal anatomy in posterior bonded restorations/provisional restorations
	OS3F	007	EF	Placement of supplemental occlusal anatomy for fine grooves
	OS3UF	007	UF	Finish supplemental occlusal anatomy in bonded/temporary restorations

By Your Side In Dentistry



BRASSLER USA
DENTAL INSTRUMENTATION