

Ultra Denture Kit





Only Available Direct from your Brasseler Representative To order call 800.841.4522 or fax 888.610.1937 In Canada call 800.363.3838 or fax 866.330.4454 Visit our website at BrasselerUSA.com



B-3736B-REVI-5022156L0-01.13-RRD

©2013 Brasseler USA. All rights reserved.

Solutions A106.11.702 Ultra Denture Kit

	Order #	Head Size	Grit	Usage
	H77E-HP	040	M	Cutter - Standard for trimming of acrylics
	H251EF-HP	060	6	Cutter - Extra-Fine for adjusting and finishing alloy/acrylic
	0634-HP	100	G	Polisher - Acrylic Extra-Oral
$ \longrightarrow $	0644-HP	100	Ø	Polisher - Acrylic Extra-Oral
	0654-HP	100	0	Polisher - Acrylic Extra-Oral

TIPS & TECHNIQUES: Cutters quickly trim acrylic and using the polishers in order will leave an ultra smooth surface.

Complimentary items: ET ProviPro (SWR22M & SWR22F)

Power Suggestion: Brasseler Forza L50K Laboratory Motor with light pressure

Sterilization: Autoclave only

Recommended Speeds: Polishers 8,000-10,000 rpms; Cutters 10,000-15,000 rpms

*This kit is also ideal for adjusting and polishing all types of acrylic appliances including retainers, splints and night guards. By Your Side In Dentistry

BRASSELER USA DENTAL INSTRUMENTATION