solutions

Drs. Ochsenbein and Becker Crown Lengthening Kit

Order Number: 232.36.500		Å	Ā	Å	Å	Å	Ā	Å	A
US Number	4	6	8	958	957	7004	7006	7008	—
Figure Number	HI	HI	HI	H207	H207	H41	H41	H41	H274
Head Size (¼ mm)	014	018	023	012	010	014	018	023	016
Head Length (mm)	_	_	_	3.0	3.0	_	_	_	3.5

Only Available Direct from your Brasseler Representative To order call 800.841.4522 or fax 888.610.1937 In Canada call 800.363.3838 or fax 866.330.4454 Visit our website at BrasselerUSA.com

By Your Side In Dentistry



B-3737B-REVI-5022159L0-01.13-RRD

Solutions[®] 232,36,500 Drs. Ochsenbein and Becker Crown Lengthening Kit

Order #	Head Size	Grit	Usage	
#4 Round Surgical	014	N/A	Removal Osteoplasty #3	
#6 Round Surgical	018	N/A	Removal Osteoplasty #2	
#8 Round Surgical	023	N/A	Removal Osteoplasty #1	
958	012	N/A	Used directly adjacent to the tooth-3mm marker - around the tooth at tooth/bone juncture	
957	010	N/A	Used directly adjacent to the tooth-3mm marker - around the tooth - interproximal	
7004	014	0	Sequentially Finishing Bone #3	
7006	018	6	Sequentially Finishing Bone #2	
7008	023	0	Sequentially Finishing Bone #1	
H274	016	0	Remove small granulated tissue tags from around the tooth	



By Your Side In Dentistry