



# BRASSELER USA®



## NEW BRASSELER USA® HAND INSTRUMENTS

### HYGIENE




SV3 Savannah-3 Scaler

- Thin contra-angle for better adaptation and enhanced posterior access
- Also recommended as instrument of choice for cement removal by leading clinician, Dr. John Kois
- Available in handles: #6, #8 (#6 shown)



SBH5/6 Barnhart 5-6 Curette

- Universal curette for medium calculus removal, also available in Rigid and Thin versions
- Available in #2, #4, #6 and #8 Cushion Grip (#6 shown)



SG11/12 Gracey 11-12 Curette

- Gracey curettes are available in a variety of designs such as Long, Mini Long, Rigid, and new LTS.
- Available in handles: #4, #6 and #8 cushion grip (shown in #6)

## DIAGNOSTIC



EX5 Explorer

- Instrument of choice for caries and calculus detection
- Available in handles: #2, #4 and #6 (shown in #6)



EX23/CP12 Expro

- Probe end also known as Marquis (3-6-9-12). Available in a variety of popular probe marking styles such as Williams (1-2-3-5-7-8-9-10) with Explorer 23.
- Available in handles: #2, #4 and #6 (shown in #6)

## RESTORATIVE - CORD PLACEMENT



CP122 - Cord Placer

- Blades turned to 0° and 90° angles from terminal shank
- Smooth and serrated versions available
- Available in handles: #2, #4, #6 and #8 cushion grip (#6 shown)



CP113 - Cord Placer

- Blades turned to a 45° angle from the terminal shank
- Smooth and serrated versions available
- Available in handles: #2, #4, #6 and #8 cushion grip (#6 shown)

## RESTORATIVE - COMPOSITE PLACEMENT



TINMINIL/R - Composite Instrument

- Mini Curved Paddle features mirror image blades 1.5 mm wide
- Titanium nitride coated to minimize composite sticking and tug back
- Available in handle #6 only



TINPF3R - Composite Instrument

- Small Paddle/Condenser features 2.0 mm wide paddle end and 1.1 mm round end condenser
- Titanium nitride coated to minimize composite sticking and tug back
- Available in handle #6 only

► Order online at [BrasselerUSA.com](http://BrasselerUSA.com) or call **800.841.4522**.  
In Canada call **800.363.3838**.

