



BRASSELER
USA®


BRiO 440

PERFORMANCE AND VALUE

The Brio440 handpiece from Brasseler USA offers a powerful yet economical alternative in highspeed handpieces. Featuring ceramic ball bearings, optic and non-optic options, 4-port water spray and swivel quick disconnect, the Brio440 provides a unique mix of performance and value in a highspeed handpiece.




- ▶ **POWERFUL 20 WATTS OF TORQUE**
- ▶ **4-PORT WATER SPRAY**
- ▶ **AVAILABLE IN OPTIC AND NON-OPTIC OPTIONS**
- ▶ **SWIVEL QUICK DISCONNECT**
- ▶ **DURABLE CERAMIC BALL BEARINGS**
- ▶ **LIGHTWEIGHT AND ERGONOMIC DESIGN**



Brio440L Optic

Item Number: 5026058U0

- Torque: 20 Watts
- Weight: 50g
- Head Height: 12.3mm
- Head Diameter: 10.5mm
- Max. Speed: 395,000 RPM
- Body: Stainless Steel



Brio440 Non-Optic

Item Number: 5026059U0

- Torque: 20 Watts
- Weight: 47g
- Head Height: 12.3mm
- Head Diameter: 10.5mm
- Max. Speed: 395,000 RPM
- Body: Stainless Steel

NL10LED

Item Number: 5021283U0

- LED Fiber Optics
- 6 Pin Tubing Connection



NL10Ti

Item Number: 5021281U0

- Fiber Optics
- 6 Pin Tubing Connection



NL11Ti

Item Number: 5021282U0

- Fiber Optics
- 5 Hole Tubing Connection



N10Ti

Item Number: 5021278U0

- Non-Optic
- 4 Hole Tubing Connection



Visit our website at BrasselerUSA.com

To order call 800.841.4522 or fax 888.610.1937. In Canada call 800.363.3838

