

NLXK™

Powerful Performance...


Up to 26 watts
of power for
smoother cutting
efficiency

Durable ceramic
ball bearings
with 3-year
warranty

Four-port
water spray for
ideal cooling

Slim neck design
for improved
visibility

LED
Optics






ComfortGRIP™ for
better feel and
protection

Lightweight
titanium with
optimum weight
distribution

360°
swivel


The NLXK high speed air handpiece from Brasseler USA is designed for maximum durability and performance. The innovative turbine design provides up to 26 watts of power, resulting in smoother cutting efficiency while reducing procedure time. Precision manufacturing and enhanced bearing technology offer improved durability backed by an industry leading three-year warranty. Equipped with ComfortGRIP coating, lightweight titanium construction and a thin neck design for improved visibility, the NLXK delivers a Powerful Performance.




NLXK™

The innovative turbine design of the NLXK provides up to **26 watts of power**, resulting in smoother cutting efficiency while reducing procedure time. Whether you are cutting through obsolete crowns or refining margins, the end result is an enhanced user experience.

Precision manufacturing and enhanced bearing technology of the NLXK lead to **improved durability, longevity and stability**. Backed by an **industry-leading three-year warranty**, the durability of the NLXK reduces your repair costs and down time while increasing your office's efficiency.



The **lightweight titanium construction** achieves optimum weight distribution while reducing hand and wrist fatigue during treatments. The **ComfortGRIP™** provides an ideal surface texture to improve grip while protecting against scratching and repeated autoclave sterilization cycles. Lastly, the **slim neck design** and **unique head angle** offer improved access and visibility.



NLXK-S - Standard Head Size

Item Number: 5026764U0

For Brasseler Coupler
 23 Watts of Power
 3-Year Warranty

NLXK-T - Torque Head Size

Item Number: 5026765U0

For Brasseler Coupler
 26 Watts of Power
 3-Year Warranty



NLXK-KS - Standard Head Size

Item Number: 5026766U0

For KaVo® MULTiflex® Coupler
 21 Watts of Power
 3-Year Warranty

NLXK-KT - Torque Head Size

Item Number: 5026767U0

For KaVo® MULTiflex® Coupler
 24 Watts of Power
 3-Year Warranty

