




A recent independent blind survey of practitioners on dental diamonds showed that Brasseler USA is rated the #1 diamond brand by a **2 to 1 margin over the closest competitor** in precision, quality and assortment and providing solutions.

INNOVATIVE DEPTH REDUCTION BURS

Developed by Dr. Robert Winter, Prosthodontist/Master Ceramist, for outcome-based preparation design.


MINIMAL
.3mm / .5mm / .7mm



RWMIN

- Preserves enamel
- 0-1 shade change
- Minor morphology changes

CONVENTIONAL
.5mm / .7mm / .9mm



RWCONV

- Preserves enamel
- 1-2 shade change
- Moderate morphology changes


EXTENSIVE
.8mm / 1.0mm / 1.2mm



RWEXT

- 2-3 shade change
- Masking of the underlying tooth structure is required
- Using an additive ceramic technique for position and morphology changes, enamel preservation is possible in the incisal half of the tooth

AVERAGE ENAMEL DEPTH



By Your Side in Dentistry



To order call 800.841.4522
or fax 888.610.1937.
Visit our website: BrasselerUSA.com



Depth reduction bur concepts and philosophy presented in Restorative Design clinics available through Spear campus learning programs: www.speareducation.com.

