

# LS SERIES Low Speed Air

**2226K1E2**  
*Durability Redefined...*

The LS Series from Brasseler USA delivers muscular performance and proven reliability for all your low speed clinical applications. You will find a unique mix of durability, strength and versatility reflected throughout this lightweight motor. The LS Series offers a variety of head and gear reduction angles, providing you with the perfect fit for your low speed procedural needs.





### LS22K Air Motor

Item Number: 5021316U0

Motor Speed Range: 0-22,000 rpm  
Connects Directly to 4 Hole Tubing  
1 Year Warranty



### LSL Head

Item Number: 5021302U0

For Latch Type Burs  
Push Button Head  
1 Year Warranty



### LSF Head

Item Number: 5021305U0

For Friction Grip Burs  
Push Button Head  
1 Year Warranty



### LS1 Sheath

Item Number: 5021296U0

1:1 Drive  
1 Year Warranty



### LS4 Sheath

Item Number: 5021298U0

4:1 Drive  
1 Year Warranty



### LS1L Sheath

Item Number: 5021300U0

1:1 Drive for Latch Type Burs  
1 Year Warranty



### LS1F Sheath

Item Number: 5021301U0

1:1 Drive for Friction Grip Burs  
1 Year Warranty



### LS1S Straight Nosecone

Item Number: 5021294U0

1:1 Drive  
For Straight Handpiece Burs and Disposable Prophy Angles  
1 Year Warranty