

ANOTHER INNOVATION FROM BRASSELER USA...

SABRECUT®

"SUPER PREMIUM" MULTI-PURPOSE CARBIDES

Alloy (Amalgam) Removal • Cavity Preparation • Oral Surgery



SMOOTH

The unique blade configuration with fine cross-cuts will generate less vibration (chatter) and leave a smoother finish.

AGGRESSIVE

The unique blade configuration and fine cross-cuts make SabreCut® carbides extremely efficient at cutting amalgam, metal, enamel and dentin.

RUGGED

Single-piece construction carbide (vs. two-piece)

- Minimizes separation
- Minimizes wear and tear on handpiece turbines.
- Reduces chair time

TIP: Two-piece conventional operative carbides are NOT recommended to remove old alloy restorations due to the high risk of separation.

END-CUTTING

With sharp, deeply ground transversing blades at the tip, SabreCut® carbides exhibit rapid end-cutting, making them ideal for restoration removal.

SabreCut® 1558SC
Compare to SS White®
Great White® No. 2.

* Internal testing performed on high carbon steel – hardness of 195 HV10 (Vickers) which is similar to the hardness of popular Nickel/Chrome/Molybdenum alloys used in crown and bridge. Complete testing available upon request.

By Your Side
in Dentistry



BRASSELER USA®
DENTAL INSTRUMENTATION




To order call 800.841.4522 or fax 888.610.1937.
In Canada call 800.363.3838 or fax 866.330.4454.
Visit our website at BrasselerUSA.com



SABRECUT®

- Multi-Purpose: Ideal for cavity preparation, alloy (amalgam) removal and oral surgery.
- Quality Manufactured: Single-piece construction carbide.
- Available in clinic packs.



Us No.	330SC	331SC*	245SC	557SC	1557SC	1558SC**	1703LSC
Head Size (1/10 mm)	008	010	008	010	010	012	021
Cutting Length (mm)	2.0	1.9	2.7	3.7	3.9	4.0	7.5
Friction Grip							
31 - FG (19 mm)	FG	FG	FG	FG	FG	FG	-
36 - FG/Surg (25 mm)	-	-	-	-	-	FG 36	-
Hall® Handpiece							
Shank #1 (44.5 mm)	-	-	-	-	-	-	HP 11
Shank #2 (51 mm)	-	-	-	-	-	-	HP 50
Shank #3 (65 mm)	-	-	-	-	-	-	HP 15
Stryker® Impaction & 100K Drills							
Shank #4 (44.5 mm)	-	-	-	-	-	-	Flat 94
Shank #5 (59 mm)	-	-	-	-	-	-	Flat 95



Scan code or visit
BrasselerUSA.com/Go/Sabrecut
 to watch our product video.

*Compare to SS White® Great White® No. 1.
 **Compare to SS White® Great White® No. 2.

By Your Side
 in Dentistry



BRASSLER USA®
 DENTAL INSTRUMENTATION



To order call 800.841.4522 or fax 888.610.1937.
 In Canada call 800.363.3838 or fax 866.330.4454.
 Visit our website at BrasselerUSA.com

