

FORZA™ Attachments

Masterpieces of Aesthetic and Structural Engineering...

Constructed of lightweight, high strength titanium, the Forza attachments are a masterpiece of aesthetic and structural engineering. The state-of-the-art gear design incorporated in the Forza attachments delivers smoother performance and increased durability. The E-Type Forza attachments are compatible with most electric motors on the market and come in a wide range of gear ratios that offer the perfect solution for your practice.





Forza F5

Item Number: 5021285U0

- 1:5 Speed Increasing
- Optic
- Speed Range: 10,000-200,000 rpm
- For Friction Grip Burs
- Four Port Spray
- 2 Year Warranty



Forza F5X

Item Number: 5021290U0

- 1:5 Speed Increasing
- Non-Optic
- Speed Range: 10,000-200,000 rpm
- For Friction Grip Burs
- Four Port Spray
- External Spray Connection
- 1 Year Warranty



Forza F1

Item Number: 5021288U0

- 1:1 Gear Ratio
- Optic
- Speed Range: 2,000-40,000 rpm
- For Latch Type Burs
- Single Port Spray
- 2 Year Warranty



Forza FIS

Item Number: 5021292U0

- 1:1 Gear Ratio
- Non-Optic
- Speed Range: 2,000-40,000 rpm
- For Straight Handpiece Burs
- Single Port Spray
- 2 Year Warranty