



# FORZA ELM™

*Changing the Way You Think About Electric Handpieces...*

The Forza ELM is arguably the most advanced clinical micromotor. The Forza ELM is a celebration of details, from the lightweight titanium motor to the LED fiber optics. With the lightest motor in its class and a powerful 3.4 Ncm of torque, this micromotor can be integrated into any dental delivery system. The Forza ELM will change the way you think about electric handpieces.

## Forza ELM Electric Micromotor

**Item Number: 5021312U0**

System Includes Motor, Cord, Control Unit and AC Adapter

Two Programmable Settings

Motor Speed Range: 2,000-40,000 rpm

Attachments Provide a System Speed Range of 100-200,000 rpm

3.4 Ncm of Torque

61 Grams Motor Weight

3 Year Warranty

LED Optics

