

# FORZA L50K™

*Uncompromising Power...*

The Forza L50K Lab Micromotor delivers uncompromising power throughout its advanced design. The cutting edge technology incorporated in the Forza L50K offers a massive 8.7 Ncm of torque and a speed range of 1,000-50,000 rpms. The lightweight, ergonomic handpiece weighs only 219 grams ensuring precise control and reduced hand fatigue.





**Forza L50K MicroMotor**  
**Item Number: 5021337U0**

Motor Speed Range: 1,000-50,000 rpm  
Maximum Torque of 8.7 Ncm  
Brushless  
Weight: 219 grams  
1 Year Warranty



**Forza L50K Complete System**  
**Item Number: 5021334U0**

System includes Motor; Cord, Control Unit,  
Foot Control and AC Adapter  
Digital Speed Display  
Easy-to-Operate Display Panel  
Foot Pedal or Hand Control Option  
Auto Cruise Function  
1 Year Warranty