

FORZA™

Clinical Electric Systems



FORZA ELS™

Muscular yet Mindful...

The Forza ELS sets the standard for clinical micromotors with a powerful 3.0 Ncm of torque in a compact, lightweight design. Smooth, quiet, and versatile, the Forza ELS features two programmable settings and an easy-to-read digital speed display. Whether cutting through obsolete crowns or refining margins, the Forza ELS delivers the power and control demanded by dentists today.

Forza ELS Electric Micromotor

Item Number: 5021306U0

System Includes Motor, Cord, Control Unit and AC Adapter

Two Programmable Settings

Motor Speed Range: 2,000-40,000 rpm

Attachments Provide a System Speed Range of 100-200,000 rpm

3.0 Ncm of Torque

76 Grams Motor Weight

3 Year Warranty





FORZA ELM™

Changing the Way You Think About Electric Handpieces...

The Forza ELM is arguably the most advanced clinical micromotor. The Forza ELM is a celebration of details, from the lightweight titanium motor to the LED fiber optics. With the lightest motor in its class and a powerful 3.4 Ncm of torque, this micromotor can be integrated into any dental delivery system. The Forza ELM will change the way you think about electric handpieces.

Forza ELM Electric Micromotor

Item Number: 5021312U0

System Includes Motor, Cord, Control Unit and AC Adapter

Two Programmable Settings

Motor Speed Range: 2,000-40,000 rpm

Attachments Provide a System Speed Range of 100-200,000 rpm

3.4 Ncm of Torque

61 Grams Motor Weight

3 Year Warranty

LED Optics



FORZA™ Attachments

Masterpieces of Aesthetic and Structural Engineering...

Constructed of lightweight, high strength titanium, the Forza attachments are a masterpiece of aesthetic and structural engineering. The state-of-the-art gear design incorporated in the Forza attachments delivers smoother performance and increased durability. The E-Type Forza attachments are compatible with most electric motors on the market and come in a wide range of gear ratios that offer the perfect solution for your practice.





Forza F5

Item Number: 5021285U0

- 1:5 Speed Increasing
- Optic
- Speed Range: 10,000-200,000 rpm
- For Friction Grip Burs
- Four Port Spray
- 2 Year Warranty



Forza F5X

Item Number: 5021290U0

- 1:5 Speed Increasing
- Non-Optic
- Speed Range: 10,000-200,000 rpm
- For Friction Grip Burs
- Four Port Spray
- External Spray Connection
- 1 Year Warranty



Forza F1

Item Number: 5021288U0

- 1:1 Gear Ratio
- Optic
- Speed Range: 2,000-40,000 rpm
- For Latch Type Burs
- Single Port Spray
- 2 Year Warranty



Forza FIS

Item Number: 5021292U0

- 1:1 Gear Ratio
- Non-Optic
- Speed Range: 2,000-40,000 rpm
- For Straight Handpiece Burs
- Single Port Spray
- 2 Year Warranty