



SPLINT SYSTEM | Developed by Brasseler USA

Designed to adjust and polish splint appliances, this system consists of E-CUTTER® instruments for trimming, diamond discs to finish and contour, a Scotch-Brite™ buff to remove tags left behind by trimming and various acrylic polishers, including ET® ProviPro™ cotton buffs to restore surface smoothness. This complete system is the ideal solution when fitting your patients with splints.

K0221

Figure No.	0674	0654	H251EF	H351EF	H79EF	365
Head Size (1/10 mm)	100	100	060	060	040	220
Cutting Length (mm)	24.0	24.0	14.0	11.0	13.0	-
Grit	Coarse	Fine	Fine	Fine	Fine	Medium
Shank Type	HP	HP	HP	HP	HP	HP
Quantity	1	1	1	1	1	1
Description	AcryliPro™ Polishing Bullet: Extra-oral polishing of acrylics		EF Cutter: Fine E-CUTTER® for trimming and finishing acrylic based appliances			Diamond Disc: Serrated for acrylics

Figure No.	SWR22M	SWR22F	RWGPP	25125HP	305
Head Size (1/10 mm)	220	220	220	250	050
Cutting Length (mm)	3.0	3.0	3.0	-	-
Grit	Medium	Fine	Medium	Fine	-
Shank Type	-	-	-	HP	HP
Quantity	1	1	1	1	3
Description	ET® ProviPro™ Polishing Buff: Polishing of acrylic and composite based appliances		Silicone Knife-Edge Polisher: Shaping and smoothing prior to final polish	Scotch-Brite™ Buff: Removal of tags created from cutters	Mounting for ET® ProviPro™ and Silicone Knife-Edge Polisher

Scotch-Brite™ is a trademark of 3M.

By Your Side
in Dentistry



BRASSLER USA
DENTAL INSTRUMENTATION



To order call 800.841.4522 or fax 888.610.1937.
In Canada call 800.363.3838 or fax 866.330.4454.
Visit our website at BrasselerUSA.com

