

SLEEP APNEA ADJUSTMENT AND POLISHING SYSTEM

Developed by Dr. Jeff Rouse



For years I have fought with adjusting and polishing the inside of night guards and sleep appliances. The soft liner that made it comfortable to wear was impossible to adjust without tearing, melting, and creating tags that were difficult to remove. The RSL Sleep Apnea Adjustment and Polishing System has resolved that problem. The choice of burs allows practitioners to adjust the appliance inside and outside and create a smooth finish even on the soft liner materials.



– Dr. Jeff Rouse

K0275

Figure No.	0674	0664	0679	0669	SWR22M	SWR22F
Head Size (1/10 mm)	100	100	055	055	220	220
Cutting Length (mm)	24.0	24.0	15.5	15.5	3.0	3.0
Grit	Coarse	Medium	Coarse	Medium	Medium	Fine
Shank Type	HP	HP	HP	HP	-	-
Quantity	1	1	1	1	1	1
Description	AcryliPro™ Polishing Bullet; Extra-oral polishing of acrylics		AcryliPro™ Polishing Point; Extra-oral polishing of acrylics		ET® ProviPro™ Polishing Buff; Polishing of acrylic and composite based appliances	




Figure No.	H79FSQ	H79FSQ	H25 IEF	25010	305
Head Size (1/10 mm)	070	040	060	100	050
Cutting Length (mm)	14.0	13.0	14.0	-	-
Shank Type	HP	HP	HP	HP	HP
Quantity	1	1	1	1	2
Description	FSQ Cutter: Fine cutter for trimming soft liner		EF Cutter: Fine E-CUTTER® for trimming and finishing acrylic based appliances	Scotch-Brite™ Buff; Smaller head size for removal of tags created by cutters	Mounting for ET® ProviPro™

Scotch-Brite™ is a trademark of 3M.

By Your Side
in Dentistry



BRASSLER USA
DENTAL INSTRUMENTATION



To order call 800.841.4522 or fax 888.610.1937.
In Canada call 800.363.3838 or fax 866.330.4454.
Visit our website at BrasselerUSA.com

