

ANOTHER INNOVATION FROM BRASSELER USA...

## HygienePro™ Ultrasonic Scaler Inserts

24 Ergonomic Scaler Inserts to Provide Greater Comfort and Efficiency for You and Your Patients!



### Wide variety of tip styles

available in both 25K and 30K frequency options to suit each user's needs

**Competitive compatibility**  
with all magnetostrictive devices  
for ease and versatility

**Patented laminate solder joint**  
to ensure outstanding durability,  
consistent spray and uniform vibration

**Ergonomic grip designs**  
feature textured non-slip surfaces for enhanced  
clinician comfort and reduced hand fatigue




#10  
UNIVERSAL

Available in slotted and direct focus water delivery options.



#100  
SLIM

Available in slotted and direct focus water delivery options.



#1000  
TRIPLE BEND

Available in slotted water delivery option only.



#3  
BEAVER TAIL

Available in slotted water delivery option only.

*Is it time to replace your ultrasonic scaler inserts?*

Did you know that only the last 1-2mm of an insert tip are the most active?

After 1mm of the tip is worn away, 25% of the tip's efficiency is lost.

After 2mm of the tip is worn away, 50% of the tip's efficiency is lost.

By Your Side  
in Dentistry



**BRASSELER USA**  
DENTAL INSTRUMENTATION



To order call 800.841.4522 or fax 888.610.1937.  
In Canada call 800.363.3838 or fax 866.330.4454.  
Visit our website at [BrasselerUSA.com](http://BrasselerUSA.com)



## HYGIENEPRO™ ULTRASONIC SCALER INSERTS - HARD GRIP



Brasseler USA® Order No.	Brasseler USA® Item No.	Hu-Friedy®	DENTSPLY® (Cavitron®)	Tip Design	Water Delivery
5023962U0	HG10U25S	UI1025K	TFI-10 25K	#10 Universal	Slotted
5023961U0	HG10U30S	UI1030K	TFI-10 30K		
5023960U0	HG10U25DF	UI25SF10	FSI-10 25K		Direct Focus
5023959U0	HG10U30DF	UI30SF10	FSI-10 30K		
5023958U0	HG100S25S	UI25K100S	TFI-100 25K	#100 Slim	Slotted
5023957U0	HG100S30S	UI30K100S	TFI-100 30K		
5023956U0	HG100S25DF	UI25SF100	FSI-100 25K		Direct Focus
5023955U0	HG100S30DF	UI30SF100	FSI-100 30K		
5023954U0	HG1000TB25S	UI100025K	TFI-1000 25K	#1000 Triple Bend	Slotted
5023953U0	HG1000TB30S	UI100030K	TFI-1000 30K		
5023963U0	HG3BT25S	UI325K	TFI-3 25K	#3 Beaver Tail	Slotted
5023964U0	HG3BT30S	UI330K	TFI-3 30K		

## HYGIENEPRO™ ULTRASONIC SCALER INSERTS - SOFT GRIP



Brasseler USA® Order No.	Brasseler USA® Item No.	Hu-Friedy®	DENTSPLY® (Cavitron®)	Tip Design	Water Delivery
5023952U0	SG10U25S	UI25SS10	TFI-10 25K	#10 Universal	Slotted
5023951U0	SG10U30S	UI30SS10	TFI-10 30K		
5023967U0	SG10U25DF	UI25SD10	FSI-10 25K		Direct Focus
5023968U0	SG10U30DF	UI30SD10	FSI-10 30K		
5023950U0	SG100S25S	UI25SS100	TFI-100 30K	#100 Slim	Slotted
5023949U0	SG100S30S	UI30SS100	TFI-100 30K		
5023948U0	SG100S25DF	UI25SD100	FSI-100 25K		Direct Focus
5023947U0	SG100S30DF	UI30SD100	FSI-100 30K		
5023965U0	SG1000TB25S	UI25SS1000	TFI-1000 25K	#1000 Triple Bend	Slotted
5023966U0	SG1000TB30S	UI30SS1000	TFI-1000 30K		
5023969U0	SG3BT25S	N/A	TFI-3 25K	#3 Beaver Tail	Slotted
5023970U0	SG3BT30S	N/A	TFI-3 30K		

## REPLACEMENT O-RINGS (PACK OF 12)



Brasseler USA Order No.	Brasseler USA Item No.	Hu-Friedy®	DENTSPLY® (Cavitron®)
5024169U0	BL12	UBLKORING	62351

Hu-Friedy® is a registered trademark of Hu-Friedy Mfg. Co., LLC.  
DENTSPLY® and Cavitron® are registered trademarks of DENTSPLY International, Inc.

By Your Side  
in Dentistry



**BRASSLER USA**  
DENTAL INSTRUMENTATION



To order call 800.841.4522 or fax 888.610.1937.  
In Canada call 800.363.3838 or fax 866.330.4454.  
Visit our website at BrasselerUSA.com

