




BRASSELER
USA®







BY YOUR SIDE

BRiO SHINE™



SINGLE-USE UNIVERSAL COMPOSITE POLISHING SYSTEM

- ▶ System adapts to all tooth surfaces
- ▶ Reduces time spent polishing
- ▶ Fewer instruments and steps


						
Order No.	BPTCOMP-1	BPTF-2	BCPCOMP-1	BCPF-2	FLBCO-1	FLBF-2
Head Size (1/10 mm)	040	040	060	060	100	100
Grit	Medium	Fine	Medium	Fine	Medium	Fine
21 - RA - Right-Angle Latch (22 mm)						
Qty/Pack	30	30	30	30	30	30

Recommended speed: 4,000-6,000 rpm, Final Polish: 2,000-3,000 rpm

- ▶ Order online at BrasselerUSA.com or call 800.841.4522.
In Canada call 800.363.3838.



©2019 Brasseler USA. All rights reserved. Printed in the U.S.A.



UNIVERSAL COMPOSITE
STEP 1

*Pre-polishing composite
restorations*



UNIVERSAL POLISHING
STEP 2

*Final, universal polisher that
provides an excellent finish
to all restorations*


B-5165-01.19-25M-PRE



BRASSELER
USA®







BY YOUR SIDE

BRIOSHINE™



SINGLE-USE UNIVERSAL CERAMIC POLISHING SYSTEM

- ▶ System adapts to all tooth surfaces
- ▶ Reduces time spent polishing
- ▶ Fewer instruments and steps

						
Order No.	BPTCER-1	BPTF-2	BCPCER-1	BCPF-2	FLBCER-1	FLBF-2
Head Size (1/10 mm)	040	040	060	060	100	100
Grit	Medium	Fine	Medium	Fine	Medium	Fine
21 - RA - Right-Angle Latch (22 mm)						
Qty/Pack	30	30	30	30	30	30

Recommended speed: 4,000-6,000 rpm, Final Polish: 2,000-3,000 rpm

- ▶ Order online at BrasselerUSA.com or call 800.841.4522.
In Canada call 800.363.3838.



©2019 Brasseler USA. All rights reserved. Printed in the U.S.A.



UNIVERSAL CERAMIC
STEP 1

Pre-polishing zirconia, lithium disilicate and porcelain/PFM restorations



UNIVERSAL POLISHING
STEP 2

Final, universal polisher that provides an excellent finish to all restorations

B-5165-01.19-25M-PRE