



BRASELER
USA®

BY YOUR SIDE

NEW!

BRIO SHINE™ STERILE



SINGLE-PATIENT POLISHING SYSTEM

- ▶ Special Feather Lite™ design adapts to all tooth surfaces
- ▶ Simple, two-step system reduces time spent polishing
- ▶ Convenient, sterile packaging

▶ Order online at BrasselerUSA.com or call 800.841.4522.
In Canada call 800.363.3838.



BRIO SHINE™ STERILE

▶ ALL CERAMIC POLISHING SYSTEM


STEP 1
Pre-Polish (Pink) for Zirconia, Lithium Disilicate, and Porcelain/PFM restorations
 Recommended Speed:
 4,000-6,000 rpm



STEP 2
Final Polish (Gray) for Zirconia, Lithium Disilicate, and Porcelain/PFM restorations
 Recommended Speed:
 2,000-3,000 rpm

Order No.	SBPTCER-1	SBPTF-2	SBCPCER-1	SBCPF-2	SFLBCER-1	SFLBF-2
Head Size (1/10 mm)	040	040	060	060	100	100
Grit	Medium	Fine	Medium	Fine	Medium	Fine
21 - RA - Right-Angle Latch (22 mm)						
Qty/Pack	25	25	25	25	25	25

▶ ALL COMPOSITE POLISHING SYSTEM


STEP 1
Pre-Polish (Green) for Composite restorations
 Recommended Speed:
 4,000-6,000 rpm



STEP 2
Final Polish (Gray) for Composite restorations
 Recommended Speed:
 2,000-3,000 rpm

Order No.	SBPTCO-1	SBPTF-2	SBCPCO-1	SBCPF-2	SFLBCO-1	SFLBF-2
Head Size (1/10 mm)	040	040	060	060	100	100
Grit	Medium	Fine	Medium	Fine	Medium	Fine
21 - RA - Right-Angle Latch (22 mm)						
Qty/Pack	25	25	25	25	25	25

▶ Order online at BrasselerUSA.com or call **800.841.4522**.
 In Canada call **800.363.3838**.

