



BRASSELER
USA®


CSERIES CARBIDES



BY YOUR SIDE

▶ OPERATIVE

- Single-piece construction for improved durability in smaller head sizes, virtually ending separation!
- Improved blade geometry in selected shapes for a faster, more aggressive cut



Order No.	¼C	½C	2C	4C	35C	330C	245C	556C	557C	169LC
Head Size (¹ / ₁₀ mm)	005	006	010	014	010	008	008	009	010	009
Cutting Length (mm)	–	–	–	–	0.9	2.0	2.8	3.3	3.9	5.0
Friction Grip										
31 - FG (19mm)	FG	FG	FG	FG	FG	FG	FG	FG	FG	FG

▶ CROWN PREP

- Available in popular chamfer and shoulder diamond bur shapes
- The aggressive blade design of the new Crown Prep Carbides efficiently removes enamel with minimum clogging



Order No.	H856CP	H856CP	H847CP	H379CP
Head Size (¹ / ₁₀ mm)	016	018	016	023
Cutting Length (mm)	8.0	8.0	8.0	4.2
Friction Grip				
31 - FG (19mm)	FG	FG	FG	FG

▶ PFM CROWN REMOVAL

- Specially designed to cut through precious and non-precious metal
- The new PFMC carbide is engineered to separate porcelain fused to metal crown in a single step



Order No.	H34C	H34C	H34LC	H34LC	PFMC
Head Size (¹ / ₁₀ mm)	010	012	012	012	012
Cutting Length (mm)	4.0	4.0	4.1	4.1	4.0
Friction Grip					
31 - FG (19mm)	FG	FG	FG	–	FG
36 - FG/Surg (25 mm)	–	–	–	FG/Surg	–

▶ TRIMMING & FINISHING

- Ideal for shaping, contouring and finishing of all composites
- Blade design reduces clogging while providing a smooth surface finish



Order No.	7801C	7802C	7803C	7404C	7406C	7408C	7404CUF	7406CUF	7408CUF	7106C	H274C
Head Size (¹ / ₁₀ mm)	009	010	012	014	018	023	014	018	023	018	016
Cutting Length (mm)	4.0	3.4	3.4	3.2	3.5	3.6	3.2	3.5	3.6	3.7	3.5

Friction Grip

31 - FG (19mm)

FG	FG	FG	FG	FG	FG	FG	FG	FG	FG	FG	FG
----	----	----	----	----	----	----	----	----	----	----	----



Order No.	7404CM	7406CM	7408CM	7901CM	7902CM	7903CM	H30011C
Head Size (¹ / ₁₀ mm)	014	018	023	009	010	012	005
Cutting Length (mm)	3.2	3.5	3.6	3.8	4.1	4.0	3.8


Friction Grip

31 - FG (19mm)

FG	FG	FG	FG	FG	FG	FG	FG
----	----	----	----	----	----	----	----

▶ SURGICAL

- Rapid cutting of bone and enamel
- Blade design ideal for efficient removal of debris



Order No.	557C	H274C	7901C	H161C	H162C	H269C
Head Size (¹ / ₁₀ mm)	010	016	009	016	016	016
Cutting Length (mm)	4.2	3.5	3.8	9.0	9.0	11.0

Friction Grip

31 - FG (19mm)


–	–	–	FG	FG	FG
---	---	---	----	----	----

36 - FG/Surg (25 mm)

FG/Surg	FG/Surg	FG/Surg	FG/Surg	FG/Surg	–
---------	---------	---------	---------	---------	---

► SPECIALTY - ENDO


- Long tapered configuration creates a funnel shape for easier access in the pulp chamber
- Non-cutting, safe-end tip widens the pulp chamber without perforation or ledging



Order No.	H269GKC	H269GKC
Head Size ($1/10$ mm)	012	016
Cutting Length (mm)	9.0	9.0
Friction Grip		
31 - FG (19mm)	FG	FG

► ORTHODONTIC DEBONDING

- Non-cutting, safe-end tip prevents damage to the gingival surface
- Designed for controlled debonding without damage to the enamel while leaving a smooth surface finish



Order No.	H283C	H48LC	H48LC	H48LCUF	7675C	H283KC
Head Size ($1/10$ mm)	012	010	012	012	016	016
Cutting Length (mm)	8.0	8.0	8.0	8.0	8.1	8.1
Friction Grip						
31 - FG (19mm)	FG	FG	FG	FG	FG	FG



Visit our website at BrasselerUSA.com

To order call 800.841.4522 or fax 888.610.1937. In Canada call 800.363.3838

