



BRASSELER
USA®

BY YOUR SIDE

Evanesce™ Flow

Flowable Restorative

► The Challenge

Today's minimally invasive dentistry demands more from flowable composites.

► The Solution

Esthetic and strong, easy to express and manipulate, Evanesce Flow is the ideal solution to minimally-invasive restorations and a vital component to your large composite cases. Evanesce Flow was engineered to complement the Evanesce Universal Restorative Composite line and features the very same chemistry.



powered by
Clinician's Choice

► Order online at BrasselerUSA.com or call **800.841.4522**.



Evanescence Flow Flowable Restorative Composite



Evanescence Flow is clearly distinguishable from composite on a radiograph for ease of diagnosis.

Whether You're Choosing for Shade, Location or Opacity—There's an Evanescence Flow For That!

- Thixotropic for excellent control, easily dispensed from the syringe, and will not slump on a vertical surface
- High strength and precise adaptation to small preparations make it ideal for superficial Class I or II restorations
- Class V restorations benefit from its no slump handling and range of opacities and shades
- Invaluable as a liner base in posterior restorations to support the dentin/adhesive seal, preventing post-op sensitivity due to its low contraction shrinkage
- Relatively low rate of shrinkage (3.6% compared with up to 5%), and high filler content (62% by weight and 48% by volume)
- Excellent shade match to the VITA Classic shade guide
- With a radiopacity of 200% AI, is clearly distinguishable from composite on a radiograph for ease of diagnosis
- Available in the most popular shades in Enamel, Universal and Dentin opacities in 1.7g syringes, Evanescence Flow will easily satisfy the majority of your flowable composite applications
- Can be used on its own or in combination with Evanescence

► Product Offering



Evanescence Flow Syringe

Shade	Opacity	Size	Order No.
A2 Dentin	90%	1.7g	323511-B
A1 Universal	85%	1.7g	323517-B
A2 Universal	85%	1.7g	323518-B
A3 Universal	85%	1.7g	323519-B
A4 Universal	85%	1.7g	323520-B
B1 Universal	85%	1.7g	323521-B
BL1 Universal	85%	1.7g	323525-B
A1 Enamel	80%	1.7g	323501-B
A2 Enamel	80%	1.7g	323502-B
A3 Enamel	80%	1.7g	323503-B
B1 Enamel	80%	1.7g	323505-B

► Order online at BrasselerUSA.com or call **800.841.4522**.

