




**BRASSELER**  
USA®

BY YOUR SIDE

## SPIRAL DIAMONDS

The new premium spiral diamonds from Brasseler USA are engineered for fast bulk reduction. The diamond-free channels create space for increased debris removal and unobstructed water cooling. The super-course spiral diamonds are available in five popular shapes and feature the perfect balance of cutting efficiency and longevity.



Order No.	TKS1	TKS2	T5856	T5856	T5379
Head Size ( $\frac{1}{10}$ mm)	012	014	016	018	023
Cutting Length (mm)	7.0	8.0	8.0	9.0	5.0
31 - FG - Friction Grip (19 mm)					
● Super-Coarse	T35005	T35005	T5856	T5856	T5379

► Order online at [BrasselerUSA.com](http://BrasselerUSA.com) or call **800.841.4522**.  
In Canada call **800.363.3838**.

