




BRASSELER
USA®

BY YOUR SIDE

VersaFlex™

Diamond-Coated Mylar Snap-On Disc & Interproximal Strip Polishing System

- Exceptional, multi-edged diamond coating provides the efficiency and versatility to polish all restorative materials including composites, ceramics, metal, and amalgam
- Color-coded to allow for easy identification of grit
- Snap-on disc feature makes changing discs easy and fast
- Flexible disc design easily conforms to the individual contours of each tooth providing unrestricted access to even hard-to-reach locations
- Corresponding mylar strips help gain access to tight interproximal areas
- Different grits on each end of strip with uncoated center for easy interproximal insertion
- For both intra-oral and extra-oral use
- Disposable, single-patient use



Clinical photographs courtesy of Douglas A. Terry, DDS.

► Order online at BrasselerUSA.com or call **800.841.4522**.
In Canada call **800.363.3838**.



VersaFlex™ Polishing System

► Product Offering


VersaFlex Snap-On Discs (100/Pack)



Head Size (1/10 mm)	12.0	10.0	12.0	10.0	12.0	10.0	12.0	10.0
Thickness (mm)	0.21		0.17		0.14		0.11	
Grit	Super-Coarse		Coarse		Medium		Fine	
Unmounted	VFD1LG	VFD1SM	VFD2LG	VFD2SM	VFD3LG	VFD3SM	VFD4LG	VFD4SM
Usage	Reduction & Removal		Smoothing & Contouring		Finishing & Pre-Polishing		High-Shine Polishing	

Recommended disc speed: 7,000 - 10,000 rpm


VersaFlex 2-Grit Interproximal Strips (100/Pack)



	Width (mm)	Length (mm)	Grit	Usage	Thickness (mm)	
VFS1W	4	150	Super-Coarse/Coarse	Removal/Smoothing	0.21	0.17
VFS2W	4	150	Medium/Fine	Pre-Polishing/High-Shine Polishing	0.14	0.11

Mandrels (5/Pk)

- Stainless steel precision mandrels for unmounted snap-on disc polishers



Length (mm)	22.0	44.5
21-RA-Right Angle Latch	310.21	—
11- HP- Straight Handpiece	—	310.11



K0315 VersaFlex Assortment Intro Pack

- 200 Snap-On Discs (25 of each size and grit)
- 50 Interproximal Strips (25 of each)
- 1 310.21 Right-Angle Latch Snap-On Mandrel
- 1 310.11 Straight Handpiece Snap-On Mandrel

► Order online at [BrasselerUSA.com](https://www.BrasselerUSA.com) or call **800.841.4522**.
In Canada call **800.363.3838**.

