

COMPATIBLE SOFTWARE

Software	Method	Special Notes
ADSTRA Dental Software	TWAIN	
Apteryx	TWAIN	
CADI V8	TWAIN	
Ciniview 10 (Instrumentarium)	TWAIN	
CareStack	TWAIN	
Carestream Imaging	TWAIN	
ClioSoft (Sota)	TWAIN	
Curve Hero	TWAIN	
DentiMax	TWAIN	
Dentrix Imaging	TWAIN	
Dexis	Custom	Requires third party TWAIN sensor module
Dolphin Imaging	TWAIN	
Dovetail	TWAIN	
Eaglesoft	TWAIN	
EVASoft Imaging Software	TWAIN	
iDentalSoft	TWAIN	
MediaDent	TWAIN	
MiPACS	TWAIN	
MOGO Cloud	Direct	
Open Dental	TWAIN	
Oryx	TWAIN	
Planmeca Romexis	TWAIN	
Progeny by Midmark	TWAIN	
Schick CDR DICOM	TWAIN	No Auto Advance
Sirona Sidexis 4	TWAIN	No Auto-Acquisition
SOPRO Imaging (Acteon)	TWAIN	
tab32	TWAIN	
TDO	TWAIN	
TigerView	TWAIN	
Gendex VixWin	TWAIN	
XDR Software	TWAIN	Version 3.2.7 or newer
XrayVision	TWAIN	

► Questions about compatibility or setup?
Contact customer service at **877.958.1625**

