

COMPATIBLE SOFTWARE		
Software	Method	Special Notes
Adstra Dental Software	TWAIN	
Apteryx™	TWAIN	
Cadi® 8	TWAIN	
Cliniview™ 10 (Instrumentarium)	TWAIN	
Carestack®	TWAIN	
Carestream™ Imaging (CS)	TWAIN	
ClioSoft (SOTA® Image)	TWAIN	
Curve Hero (Curve Dental®)	TWAIN	
Dentimax™	TWAIN	
Dentrix™ Imaging	TWAIN	
Dexis™	Custom	Requires third party TWAIN sensor module
Dolphin ImagingPlus™	TWAIN	
Dovetail	TWAIN	
Eaglesoft™	TWAIN	
EvaSoft Imaging Software	TWAIN	
EZ Dent-I Imaging (Vatech)	TWAIN	
iDentalSoft™	TWAIN	
Mediadent™	TWAIN	
MiPACS™	Direct	
Open Dental™	TWAIN	
Oryx™	TWAIN	
Planmeca Romexis™	TWAIN	
Progeny™ by Midmark	TWAIN	
Schick CDR Dicom	TWAIN	No Auto Advance
Sirona SIDEXIS™ 4	TWAIN	No Auto-Acquisition
SoPro Imaging™ (Acteon)	TWAIN	
Tab32™	TWAIN	
TD0™	TWAIN	
TigerView™	TWAIN	
VixWin™ (Gendex)	TWAIN	
XDR™ Software	TWAIN	Version 3.2.7 or newer
XrayVision™	TWAIN	
Mogo Cloud™	TWAIN	

Adstra Dental Software is a registered trademark of Adstra, Inc. Apteryx™ and XrayVision™ are registered trademarks of Apteryx, LLC. Cadi® 8 is a registered trademark of Synca Marketing Inc. Cliniview™ 10 (Instrumentarium) is a registered trademark of Palodex Group Oy OY FINLAND. Carestack® is a registered trademark of Good Methods Global Inc. Carestream™ Imaging (CS) is a registered trademark of Carestream Health, Inc. ClioSoft (SOTA® Image) is a registered trademark of Precision Optics, Inc. Curve Hero (Curve Dental®) is a registered trademark of Curve Dental Ltd. Dentimax™ is a registered trademark of HENRY SCHEIN MEDICAL SYSTEMS, INC. Dentrix™ Imaging is a registered trademark of HS™, LLC. Dexis™ is a registered trademark of Dental Imaging Technologies Corporation. Dolphin ImagingPlus™ is a registered trademark of Patterson Dental Holdings, Inc. Dovetail is a registered trademark of DOVETAIL DENTAL ASSOCIATES LIMITED LIABILITY COMPANY. Eaglesoft™ is a registered trademark of Patterson Companies, Inc. EvaSoft Imaging Software is a registered trademark of SMK Imaging LLC. EZ Dent-I Imaging (Vatech) is a registered trademark of VATECH CO., LTD CORPORATION REPUBLIC OF KOREA. DentalSoft™ is a registered trademark of iDentalSoft, Inc, a Delaware corporation. Mediadent™ is a registered trademark of MultiMedia Dental Systems, Inc. MiPACS™ is a registered trademark of Lead Technologies, Inc. CORPORATION NORTH CAROLINA. Open Dental™ is a registered trademark of Open Dental Software. Oryx™ is a registered trademark of Oryx Dental Software. Planmeca Romexis™ is a registered trademark of Planmeca Oy CORPORATION FINLAND. Progeny™ by Midmark is a registered trademark of Midmark Corporation. Schick CDR Dicom and Sirona SIDEXIS™ 4 are registered trademarks of Dentsply Sirona. SoPro Imaging™ (Acteon) is a registered trademark of Acteon Group. Tab32™ is a registered trademark of Integrated Charts, Inc. TD0™ is a registered trademark of MedMark, LLC. TigerView™ is a registered trademark of Televere Systems, LLC LIMITED LIABILITY COMPANY. VixWin™ (Gendex) is a registered trademark of Dental Imaging Technologies Corporation. XDR™ Software is a registered trademark of Cyber Medical Imaging, Inc.. Mogo Cloud™ is a registered trademark of Mogo, Inc..

► Questions about compatibility or setup?
Contact customer service at **877.958.1625**



BRASSELER
USA®