



BRASSELER
USA®

CSERIES CARBIDES

BY YOUR SIDE



▶ OPERATIVE

- Single-piece construction for improved durability in smaller head sizes, virtually ending separation
- Improved blade geometry in selected shapes for a faster, more aggressive cut

										
Order No.	1/4C	1/2C	2C	4C	35C	330C	245C	556C	557C	169LC
Head Size ($\frac{1}{10}$ mm)	005	006	010	014	010	008	008	009	010	009
Cutting Length (mm)	—	—	—	—	0.9	2.0	2.8	3.3	3.9	5.0
Friction Grip										
31 - FG (19mm)	FG	FG	FG	FG	FG	FG	FG	FG	FG	FG

▶ CROWN PREP

- Available in popular chamfer and shoulder diamond bur shapes
- Aggressive blade design efficiently removes enamel with minimum clogging

				
Order No.	H856CP	H856CP	H847CP	H379CP
Head Size ($\frac{1}{10}$ mm)	016	018	016	023
Cutting Length (mm)	8.0	8.0	8.0	4.2
Friction Grip				
31 - FG (19mm)	FG	FG	FG	FG

▶ PFM CROWN REMOVAL

- H34 series is specially designed to cut through precious and non-precious metal
- PFMC carbide is engineered to separate porcelain fused to metal crown in a single step

					
Order No.	H34C	H34C	H34LC	H34LC	PFMC
Head Size ($\frac{1}{10}$ mm)	010	012	012	012	012
Cutting Length (mm)	4.0	4.0	4.1	4.1	4.0
Friction Grip					
31 - FG (19mm)	FG	FG	FG	—	FG
36 - FG/Surg (25 mm)	—	—	—	FG/Surg	—

▶ TRIMMING & FINISHING

- For shaping, contouring and finishing of all composites
- Blade design reduces clogging while providing a smooth surface finish

											
Order No.	7801C	7802C	7803C	7404C	7406C	7408C	7404CUF	7406CUF	7408CUF	7106C	H274C
Head Size ($\frac{1}{10}$ mm)	009	010	012	014	018	023	014	018	023	018	016
Cutting Length (mm)	4.0	3.4	3.4	3.2	3.5	3.6	3.2	3.5	3.6	3.7	3.5
Friction Grip											
31 - FG (19mm)	FG	FG	FG	FG	FG	FG	FG	FG	FG	FG	FG

							
Order No.	7404CM	7406CM	7408CM	7901CM	7902CM	7903CM	H30011C
Head Size ($\frac{1}{10}$ mm)	014	018	023	009	010	012	005
Cutting Length (mm)	3.2	3.5	3.6	3.8	4.1	4.0	3.8
Friction Grip							
31 - FG (19mm)	FG	FG	FG	FG	FG	FG	FG


▶ SURGICAL

- Rapid cutting of bone and enamel
- Blade design efficiently removes debris

							
Order No.	557C	H274C	7901C	H161C	H162C	H269C	
Head Size ($\frac{1}{10}$ mm)	010	016	009	016	016	016	
Cutting Length (mm)	4.2	3.5	3.8	9.0	9.0	11.0	
Friction Grip							
31 - FG (19mm)	—	—	—	FG	FG	FG	
36 - FG/Surg (25 mm)	FG/Surg	FG/Surg	FG/Surg	FG/Surg	FG/Surg	—	

SPECIALTY - ENDO


- Long tapered configuration creates a funnel shape for easier access in the pulp chamber
- Non-cutting, safe-end tip widens the pulp chamber without perforation or ledging



Order No.	H269GKC	H269GKC
Head Size ($\frac{1}{10}$ mm)	012	016
Cutting Length (mm)	9.0	9.0
Friction Grip		
31 - FG (19mm)	FG	FG

ORTHODONTIC DEBONDING


- Non-cutting, safe-end tip prevents damage to the gingival surface
- Designed for controlled debonding without damage to the enamel while leaving a smooth surface finish



Order No.	H283C	H48LC	H48LC	H48LCUF	7675C	H283KC
Head Size ($\frac{1}{10}$ mm)	012	010	012	012	016	016
Cutting Length (mm)	8.0	8.0	8.0	8.0	8.1	8.1
Friction Grip						
31 - FG (19mm)	FG	FG	FG	FG	FG	FG

PERI-IMPLANTITIS

- Developed with Dr. Dennis Tarnow and Dr. Paul Rosen for treatment of peri-implantitis
- Features an extra-long, 30mm surgical-length, FG shank



Order No.	TR1	TR1UF	TR2	TR2UF	TR3	TR3UF	TR4	TR4UF	TR5	TR5UF
Head Size ($\frac{1}{10}$ mm)	014	014	023	023	014	014	023	023	016	016
Cutting Length (mm)	8.0	8.0	8.0	8.0	3.1	3.1	4.1	4.1	3.5	3.5
Friction Grip										
30 - XL FG Surg (30mm)	XL FG Surg	XL FG Surg	XL FG Surg	XL FG Surg	XL FG Surg	XL FG Surg	XL FG Surg	XL FG Surg	XL FG Surg	XL FG Surg

Visit our website at BrasselerUSA.com

To order call 800.841.4522 or fax 888.610.1937. In Canada call 800.363.3838

